RNA Isolation Protocol – Damla D. Bilgin-IGB, UIUC.

Materials:

Acid Phenol:Chloroform: Ambion Cat# 9722 Ammonium Thiocyanate: Sigma A7149 Chloroform (HPLC Grade): Fisher Cat# C606SK-1 Glycerol: Fisher G33-500 Guanidine Thiocyanate: Ambion Cat# 9422 Phenol: Sigma P4557 3M Potassium Acetate pH 5.5: Ambion Cat# 9610 3M Sodium Acetate pH 5.5: Ambion Cat# 9740

HMT Preparation Protocol:

Added		Conc.	20mL	40mL	100mL	200mL
	Phenol Solution	38%	7.6mL	15.2mL	38.0mL	76.0mL
	Guanidine Thiocyanate	0.8M	1.89g	3.78g	9.45g	18.9g
	Ammonium Thiocyanate	0.4M	0.638g	1.276g	3.19g	6.38g
	3M Sodium Acetate pH 5.5	0.1M	0.668mL	1.336mL	3.34mL	6.68mL
	Glycerol	5%	1.0mL	2.0mL	5.0mL	10.0mL

Add nuclease-free water to complete to volume.

Protocol

- Pre-spin empty 50¹mL Phase Lock Gel-Heavy tubes (PLG) briefly to collect gel on tube bottoms (1500g for 1 min @ 4°C is sufficient to collect gel at tube bottoms).
- 2. Pre-cool mortar with liquid nitrogen. Add excess liquid nitrogen so it does not fully evaporate.
- 3. Add frozen tissue to mortar containing liquid nitrogen.
- 4. Grind to a fine powder, carefully adding more nitrogen, if necessary, as it evaporates to ensure tissue does not thaw². *Be careful not to splash tissue out of the mortar*
- 5. Once ground sufficiently, let nitrogen evaporate off and add 6¹mL of HMT and keep on grinding.
- 6. Cut the end off a pipette tip with a clean razor and transfer the sample to 50¹mL PLG tube.
- 7. Add $2^{1,3}$ mL HMT to the mortar again and rinse the mortar and transfer the liquid to the same PLG tube. (Total of 8^4 mL)
- 8. Add 1/5 volume of chloroform (1/5 * Total vol. of HMT added)
- 9. Mix by rocking⁶ and incubate at room temperature for 5 minutes.
- 10. Centrifuge swing out rotor program 2 (2000g for 10min @ 4°C)

- Transfer and measure aqueous phase to a new 50¹mL PLG tube add cold 3M Potassium Acetate pH 5.5 (1mL Potassium Acetate for every 3mL of aqueous layer)
- 12. Mix by swirling⁶
- 13. Add equal volume of acid phenol:chloroform (1mL / 1mL aqueous layer)
- 14. Mix by swirling⁶
- 15. Incubate on ice for 1hr
- 16. Centrifuge swing out rotor, program 2 (2000g for 10min @ 4°C)
- 17. Transfer⁷ and measure⁵ the aqueous phase to 'oak ridge' tube
- 18. Add 0.25 mL of 0.8 M sodium citrate/1.2 M NaCl per 1 mL of aqueous phase transferred.
- 19. Mix by swirling⁶
- 20. Add 0.25mL of absolute EtOH / 1mL aqueous (0.25*aqueous phase transferred)
- 21. Mix by gentle inversion⁶ and incubate at -20 °C for at least 30 minutes (overnight incubation produces higher yields).
- 22. Centrifuge fixed rotor, program 3 (11,000g for 10min @ 4°C)
- 23. Decant the supernatant (Being careful not to lose pellet)
- 24. Add 2mL of 70% EtOH, vortex
- 25. Centrifuge, program 3 (11,000g for 10min @ 4°C)
- 26. Discard supernatant (being careful not to lose pellet)
- 27. Repeat 70% EtOH wash
- 28. Centrifuge, program 3 (11,000g for 10min @ 4°C)
- 29. Discard supernatant (being careful not to lose pellet)
- 30. Air dry the pellet
- 31. Dissolve in RNase free water (200 µL)

Notes:

- 1-depends on starting material
- 2 dark green = thawed
- 3 Add appropriate amount of HMT until glue-like consistency
- 4 Based on amount added of HMT added in step 7
- 5 Volume should be equal to amount HMT added
- 6- DO NOT SHAKE!

7 – If Phase Lock broke or there seems to be gel in sample then transfer sample to a new PLG tube and add equal amounts of acid phenol:chloroform and centrifuge again with a water a counter balance.

Name: Ouick Glance Flow Chart:		Da		
Quick Glunce Flow Churt.	Sample	Sample	Sample	Sample
-Add HMT to mortat to transfer and wash				
-Add 1/5 volume chloroform				
-Incubate 5 min while rocking-Centrifuge Program 2-Transfer and measure aqueous to new PLG				
-Add 1/3 volume KAc and swirl				
-Add 1 vol. of acid phenol:chloroform - sw	irl			
-Incubate 1 hour on ice-Centrifuge Program 2-Transfer and measure aqueous to oak ridge				
-Add ¹ / ₄ volume of NaCitrate/NaCl and swirl	l			
-Add ¹ / ₄ volume of 100% EtOH				
-Gently Invert and incubate at -20°C overnig -Centrifuge Program 3	ht			
-Add 2mL 70% EtOH -Centrifuge Program 4		2n	nL	
-Discard supernatant -Add 2mL 70% EtOH -Centrifuge Program 4 -Discard Supernatant		2n	nL	
-Air Dry pellet -Dissolve in 200µL of nuclease-free water		200)μL	
Sample # \rightarrow	Sample #	_>		
Sample # \rightarrow	Sample #	_→		